

VERSION 1 : 2023

Prescription Drug Monitoring Reference Guide

CompuNet Tox Line:
1-833-TOX-4CCL (1-833-869-4225)

Monday - Friday
8:00 AM - 4:30 PM

 **CompuNet**[®]
CLINICAL LABORATORIES

Prescription Drug Monitoring Reference Guide

						Source of Drug Metabolites		
DRUG/DRUG CLASS	SCREEN METHOD OF ANALYSIS	INITIAL SCREEN CUTOFF (ng/ml)	LCMS CONFIRMATIONS	DRUG QUANTITATIONS/ CONFIRMATIONS	ADDITIONAL METABOLITES/ COMPOUNDS OFTEN PRESENT	GENERIC NAME	COMMON TRADE/ STREET NAME	URINE DETECTION TIME
Alcohol Metabolites			500	Ethyl Glucuronide (EtG)	Ethyl Sulfate (EtS)			Up to 80 hours
			100	Ethyl Sulfate (EtS)	Ethyl Sulfate (EtS)			
Amphetamines	EIA	500	200	Amphetamine	None	Amphetamine	Adderall®, Dexedrine®, Dextrostat®, Vyvanse®	1-3 Days
			200	Methamphetamine	Amphetamine	Methamphetamine	Desoxyn®	
			200	MDMA	MDA	MDMA	Ecstasy	
			200	MDA	None	MDA	Ecstasy/Sally	
			200	MDEA	MDA	MDEA	Ecstasy	
Amphetamines with Methamphetamine d/l Isomers	EIA			Amphetamine	None	Amphetamine	Adderall®, Dexedrine®, Dextrostat®, Vyvanse®	1-3 Days
			200	Methamphetamine	Amphetamine	Methamphetamine	Desoxyn®	
				d-Methamphetamine	None	d-Methamphetamine	Desoxyn®	
				l-Methamphetamine	None	l-Methamphetamine	Nasal Inhaler®, Selegiline®	
Barbiturates	EIA	200	100	Amobarbital	(Secobarbital Combination)	Amobarbital	Amytal®, Tuinal® - (Amobarbital/ Secobarbital Combination)	1-3 Days Short Acting; > 3 Long Acting
			100	Butalbital	None	Butalbital	Esgic®, Fioricet®, Fiorinal®, Sandoptal®	
			100	Pentobarbital	None	Pentobarbital	Nembutal®	
			100	Secobarbital	(Amobarbital Combination)	Secobarbital	Seconal®, Tuinal® - (Amobarbital/ Secobarbital Combination)	
			100	Phenobarbital	None	Phenobarbital	Luminal®	

DRUG/DRUG CLASS	SCREEN METHOD OF ANALYSIS	INITIAL SCREEN CUTOFF (ng/ml)	LCMS CONFIRMATIONS	DRUG QUANTITATIONS/ CONFIRMATIONS	ADDITIONAL METABOLITES/ COMPOUNDS OFTEN PRESENT	Source of Drug Metabolites		URINE DETECTION TIME
						GENERIC NAME	COMMON TRADE/ STREET NAME	
Benzodiazepine	EIA	100	25	Alphahydroxyalprazolam	None	Alprazolam	Xanax®	1-3 Days Short Acting; > 3 Long Acting
			50	Alphahydroxymidazolam	None	Midazolam	Versed®	
			50	Alphahydroxytriazolam	None	Triazolam	Halcion®	
			25	Aminoclonazepam	None	Clonazepam	Klonopin®	
			50	Hydroxyethylflurazepam	None	Flurazepam	Dalmane®	
			50	Lorazepam	None	Lorazepam	Ativan®	
			50	Nordiazepam	Oxazepam, Temazepam	Chlordiazepoxide, Clorazepate, Diazepam	Chlordiazepoxide-Librium®, Clorazepate-Tranxene®, Diazepam-Valium®	
			50	Oxazepam	Nordiazepam, Temazepam	Chlordiazepoxide, Clorazepate, Diazepam, Oxazepam, Temazepam	Chlordiazepoxide-Librium®, Clorazepate-Tranxene®, Diazepam-Valium®, Oxazepam-Librax®, Serax®, Temazepam-Restoril	
		50	Temazepam	Oxazepam, Nordiazepam	Diazepam, Temazepam	Restoril®, Valium®		
Buprenorphine	EIA	5	5	Buprenorphine	Norbuprenorphine	Buprenorphine	Suboxone®, Zubsolv®, Bunavil®	1-3 Days
			5	Norbuprenorphine	None	Norbuprenorphine	Subutex®, Bunavil®, Buprenex®, Butrans®	
			5	Naloxone	None	Naloxone	Suboxone®, Bunavil®, Zubsolv®	
Cocaine Metabolite	EIA	150	100	Benzoyllecgonine	None	Cocaine	Coke	1-3 Days
Fentanyl	EIA	1	0.5	Fentanyl	Norfentanyl	Fentanyl	Abstral®, Actiq®, Duragesic®, Fentora®, Ionsys®	1-2 Days
			1	Norfentanyl	None	Fentanyl		
Gabapentin			1000	Gabapentin	None	Gabapentin	Fantrex, Gabarone, Neurontin®	1-3 Days
Heroin Metabolite			10	6-Acetylmorphine (6-AM)	None	Heroin	Heroin	1 Day
Marijuana Metabolite	EIA	20	10	Marijuana Metabolite	None	Marijuana	Marinol®, Dronabinol®	1-3 Days Casual Use; > 3 Chronic Use
Methadone Metabolites (EDDP)	EIA	150	100	MAD	EDDP	Methadone	Dolophine	1-3 Days
		150	100	EDDP	None			
Methylphenidate			100	Ritalinic Acid	None	Methylphenidate	Ritalin®, Concerta®	1-3 Days
Naltrexone	EIA		5	Naltrexone	6-beta-naltrexol	naltrexone	ReVia, Vivitrol	1-3 Days

DRUG/DRUG CLASS	SCREEN METHOD OF ANALYSIS	INITIAL SCREEN CUTOFF (ng/ml)	LCMS CONFIRMATIONS	DRUG QUANTITATIONS/ CONFIRMATIONS	ADDITIONAL METABOLITES/ COMPOUNDS OFTEN PRESENT	Source of Drug Metabolites		URINE DETECTION TIME
						GENERIC NAME	COMMON TRADE/ STREET NAME	
Nicotine & Metabolites	EIA	500	15					
Opiates	EIA	150	50	Codeine	Morphine (Hydrocodone-minor metabolite)	Codeine	Tylenol #3, Tylenol #4, Tylenol with Codeine, Acetaminophren with Codeine, Phenaphen with Codeine, Phenergen with Codeine, Phenameth with Codeine, Promethasine with Codeine	1-3 Days
		150	50	Morphine	(Hydromorphone-minor metabolite)	Morphine	Avinza, Kadian, MS Contin, Roxanol	
		150	50	Hydrocodone	Hydromorphone, Norhydrocodone	Hydrocodone	Hycodan, Lorcet, Lortab, Norco, Vicodin	
		150	50	Hydromorphone	None	Hydromorphone	Dilaudid	
		150	50	Norhydrocodone	Hydromorphone	Hydrocodone	Hycodan, Lorcet, Lortab, Norco, Vicodin	
Oxycodone and Metabolites	EIA	100	50	Noroxycodone	Oxycodone, Oxymorphone	Oxycodone	Endocet, Magnacet, OxyCotiin, OxyFast, OxyIR, Percocet, Percodan, Perolox, Roxicodone, Tylox	1-3 Days
		100	50	Oxycodone	Nnoroxycodone, Oxymorphone	Oxycodone	Endocet, Magnacet, OxyContin, OxyFast, OxyIR, Percocet, Percodan, Perolox, Roxicodone, Tylox	
		150	50	Oxymorphone	None	Oxymorphone	Opana, Numorphan	
Phencyclidine	EIA	25	5	Phencyclidine	None	Phencyclidine	PCP	1-3 Days
Pregabalin			1000	Pregabalin	None	Pregabalin	Lyrica®	1-3 Days
Propoxyphene	EIA	300	200	Norpropoxyphene	None	Propoxyphene	Darvon®, Davocet®, Dolene®, Progesic®, Trucet®	1-3 Days
Tapentadol			50	Tapentadol	Nortapentadol	Tapentadol	Nucynta®	1-3 Days
			50	Nortapentadol	None	Tapentadol	Nucynta®	
Tramadol			100	Tramadol	Desmethyltramadol	Tramadol	Rybix®, Ryzolt®, Ultram®	1-3 Days
			100	Desmethyltramadol	None	Tramadol	Rybix®, Ryzolt®, Ultram®	
Tricyclic Antidepressants	EIA	300	25	Amitriptyline	Nortriptyline	Amitriptyline, Nortriptyline	Elavil®, Endep®	1-3 Days
	EIA	300	25	Nortriptyline	None	Nortriptyline	Aventyl, Pamelor	
Zolpidem			20	Zolpidem	None	Zolpidem	Ambien®, Edluar®, Intermezzo®, Zolpimist®	1-3 Days
Zolpidem Metabolites			20					

ABOUT COMPUNET CLINICAL LABORATORIES

CompuNet has provided clinical and pathology lab services throughout Southwest Ohio since 1986. Our clients range from a health system with a Level 1 Trauma center, regional hospitals and a large physician network to smaller, one provider independent offices. CompuNet has a large patient service center network throughout the region and accepts all healthcare provider orders. CompuNet also offers direct access to patients for some testing. We offer over 2,000 tests ranging from glucose testing to complex cancer diagnostics.

Our mission is to provide quality diagnostic laboratory services to improve health outcomes within the communities we serve.

**LAB RESULTS INSPIRE
A HEALTHIER YOU™**

compunetlab.com



For assistance with interpreting Prescription Drug Testing results, please contact a CompuNet Clinical Laboratories Specialist at 1-833-TOX-4CCL or 1-833-869-4225 M-F 8:00 AM - 4:30 PM

